

STRONGARM® healthcare

OneCart™

A flexible, mobile solution for bedside patient care, telehealth applications, and more.



UV Sterilization Approved

Limited Lifetime Warranty

- **Modular, Customizable Components**
- **Ergonomic and Adjustable Design**
- **Concealed Cable Management**


- **Industrial-Grade for Stability**
- **Medical-Grade Locking Casters**
- **Antimicrobial Finish**

OneCart™




Proudly built in the U.S.A.

Modular and Ergonomic Design
Industrial Grade
UV Sterilization Approved
Limited Lifetime Warranty




MightyMount Features:

- Repositions effortlessly
 - 100% counter-balanced
 - Secure lock-in-place lever
 - 12" of vertical articulation
 - Dual-axis pivot action
 - Tilting display mount
 - Standard VESA bracket
- See more on our website...




DISPLAYS AND KEYBOARD ENCLOSURES/TRAYS MAY ALSO BE MOUNTED DIRECTLY TO TRACK

OneCart. Unlimited possibilities...



OneCart's modular and customizable design makes it the ideal mobile solution for bedside patient care and telehealth applications.


OneCart™ Modular Components & Accessories



**MIGHTYMOUNT
ARTICULATING ARM
WITH KEYBOARD
ENCLOSURE, CABLE
MANAGEMENT, AND
VESA DISPLAY MOUNT**

MightyMount Features:

- Repositions effortlessly
- 100% counter-balanced
- Secure lock-in-place lever
- 12" of vertical articulation
- Dual-axis pivot action
- Tilting display mount
- Standard VESA bracket



**KEYBOARD ENCLOSURE WITH WORK SURFACE AND
INTEGRATED USB HUB - AVAILABLE IN 18" or 22" WIDTH
(OPTIONAL LEFT/RIGHT SLIDING MOUSE TRAY)**



**VESA DISPLAY MOUNT
WITH TILT**



**VESA DISPLAY MOUNT
WITH TILT & SWIVEL**



**WEBCAM MOUNT
WITH TILT**



PTZ CAMERA MOUNT



**SPEAKER / AUDIO
CONTROL MODULE MOUNT**




**MULTI-RECEPTABLE
120V JUNCTION BOX**



POWER CORD CLIP



LOCKABLE UTILITY BOX




ACCESSORY TRAY



**UNINTERRUPTIBLE
POWER SUPPLY (UPS)**



**OUR INNOVATIVE CAM LOCK
SYSTEM ALLOWS MOUNTS AND
ACCESSORIES TO BE MOUNTED AND
REPOSITIONED INDEPENDENTLY,
WITHOUT HAVING TO REMOVE
OTHER ITEMS FROM THE TRACK!**



**ONECART'S VERSATILE DUAL-TRACK SYSTEM
ALLOWS FRONT AND REAR MOUNTING**

**Need a customized mount or accessory for your OneCart solution?
Contact our sales team for details!**

Strongarm Healthcare
425 Caredean Drive
Horsham, PA 19044

Toll-free: (800) 778-7901
Phone: (215) 443-3400
Fax: (215) 443-3002

