

Supporting Excellence in Patient Care
Delivery with Care Intelligence™

STRONGARM®
+ healthcare

*Innovative, reliable solutions that enhance efficiency, improve ergonomics,
and empower healthcare professionals to deliver exceptional patient care*



About Strongarm Healthcare

Design-First Approach

Our design-first approach ensures every solution integrates seamlessly into your workflow. Whether you're planning additional patient rooms, upgrading existing spaces, or standardizing equipment across an entire health system, we provide tailored solutions that improve efficiency, enhance ergonomics, and elevate patient care.

Engineering Expertise

Backed by more than 30 years of engineering expertise, we design ergonomic workstation solutions and mounting systems that deliver reliable performance, exceptional durability, and lasting value in demanding healthcare environments.

In-House Manufacturing

Every Strongarm Healthcare product is engineered, manufactured, and assembled at our Philadelphia, Pennsylvania facility, ensuring exceptional quality, dependable lead times, and reliable performance. Our vertically integrated manufacturing process enables rapid customization, faster delivery, and industry-leading support for every product we build.

Trusted U.S. Healthcare Partner

Since 1990, Strongarm has proudly designed and manufactured innovative healthcare solutions in the USA from our 100,000-square-foot Philadelphia, Pennsylvania facility. More than a hardware supplier, we're a trusted partner, collaborating with you from concept through installation to deliver seamless equipment integration and an exceptional clinical experience.



Care Intelligence™ is the thoughtful attention to every detail, inspired by a genuine understanding of the people who rely on our products.

Healthcare professionals bring extraordinary care and precision to everything they do. At Strongarm Healthcare, we apply that same commitment to everything we create.

From design and engineering to manufacturing, assembly, and support, Care Intelligence guides every step of our process. We thoughtfully design and carefully craft every product to deliver exceptional quality, lasting reliability, and dependable performance.

Because our commitment to care extends well beyond manufacturing, we back our products with a Lifetime Warranty.

Care Intelligence. Thoughtfully designed. Carefully built. Made to last.

Product Line-Up

STRONGARM[®]
+healthcare



CleanMount TRACK MOUNTED WALL ARM

CleanMount features a fully enclosed, pre-cabled passageway that eliminates external cable clutter, while smooth, non-porous surfaces make cleaning quick and effective. Multiple pivot points and counterbalanced vertical articulation allow the display and keyboard to be effortlessly repositioned for optimal ergonomics and user comfort. Constructed from durable cast aluminum, CleanMount is built to withstand demanding healthcare environments and is ideal for hospitals, clean rooms, and point-of-care settings where cleanliness, durability, and reliability are essential.




MightyMount TRACK MOUNTED WALL ARM

MightyMount combines ergonomic performance with exceptional durability. Dual-axis pivot points and smooth, counterbalanced vertical articulation allow the display and keyboard to be effortlessly repositioned for optimal comfort and ergonomics. Its antimicrobial, easy-to-clean hygienic finish helps support infection control, while a wide range of customization options allows each wall arm system to be tailored to your specific point-of-care application. Constructed from industrial-grade die-cast aluminum, MightyMount delivers lasting reliability in even the most demanding healthcare environments.

CleanMount

Customization & Accessories




Enhance your CleanMount system with a wide range of customizable accessories, including PC holders and enclosures, horizontal arm extensions, dual-display mounts, barcode scanner brackets, RFID scanner holders, card readers, webcams, and a variety of keyboard enclosures and trays. Designed for seamless integration, these options allow you to configure each system to meet your specific workflow and application requirements.

DISPLAY, RFID READER, SCANNER, KEYBOARD, AND MOUSE NOT INCLUDED.

*OPTIONAL ACCESSORY CHOICES

MightyMount

Customization & Accessories



Enhance your MightyMount system with a wide range of customizable accessories, including PC holders and enclosures, horizontal arm extensions, dual-display mounts, barcode scanner brackets, RFID scanner holders, card readers, webcams, and a variety of keyboard enclosures and trays. Designed for seamless integration, these options allow you to configure each system to meet your specific workflow and application requirements.


DISPLAY, RFID READER, SCANNER, KEYBOARD, AND MOUSE NOT INCLUDED.

*OPTIONAL ACCESSORY CHOICES

OneCart MOBILE & MODULAR WORKSTATION


OneCart is Strongarm Healthcare's flexible, modular, and fully customizable mobile workstation, designed for bedside patient care, telehealth, and a wide range of clinical applications. Backed by the industry's only limited lifetime warranty, OneCart delivers exceptional durability, adaptability, and long-term value.

OneCart. Unlimited Possibilities.



Telehealth Stations MODULAR & CUSTOMIZABLE

Our modular telehealth solutions are customized to meet your application requirements, with a wide range of configurable options to support virtually any clinical workflow. Available options include VESA display mounts, webcams, PTZ cameras, speakers, audio control modules, keyboard enclosures and trays, PCs, power supplies, barcode scanners, card readers, and more. Every telehealth solution is UV sterilization compatible and backed by the industry's only limited lifetime warranty.






Medical Device Mountings HIGHLY CUSTOMIZABLE

Strongarm Healthcare designs and manufactures mounting solutions for a wide range of clinical medical devices. Available mounting configurations include wall mount, ceiling mount, countertop mount, and track or rail mount to support virtually any care environment. With rapid engineering, in-house manufacturing, and fast delivery, we're the trusted partner for your next fit-out, point-of-care upgrade, or custom device mounting solution.




FlexMount AV and **FlexMount AV Max** INDUSTRIAL-GRADE AV MOUNTINGS

Engineered for demanding clinical and AV environments, FlexMount AV combines industrial-grade durability with ergonomic performance. VESA compatible and designed to support a wide range of display sizes, it delivers secure installation, smooth adjustability, and reliable positioning for enhanced user comfort and accessibility.



FlexMount AV

FlexMount AV Max




Customized Solutions

EXPERIENCE + ENGINEERING EXCELLENCE

As a pioneer in healthcare workstation and mounting solutions, Strongarm Healthcare brings more than three decades of engineering expertise to every project. Since 1990, we've designed, manufactured, and delivered more than **580,000 products**, including **410,000+ unique designs**, each tailored to our customers' unique applications. Backed by extensive in-house engineering and manufacturing capabilities, we create purpose-built solutions that integrate seamlessly into your workflow and support long-term performance.





CleanMount™




clean (klēn)

1. adjective

Something that is clean is free from dirt or contamination.

Synonyms: hygienic, sterile





Design Flow Chart



Client Assessment

Every project begins with a thorough assessment of your needs, specifications, workflow, and environment, ensuring a solution tailored to your application.

Collaborative Approach




Our engineering team works closely with healthcare planners, architects, and stakeholders throughout each project, providing transparency, flexibility, and alignment with your vision from concept to installation.

Advanced Design Resources

Using advanced CAD modeling and simulation tools, our engineers develop precise, high-performance designs with speed and accuracy.

Prototyping & Testing

From rapid prototyping to material testing and compliance verification, we validate every design to ensure lasting performance, durability, and reliability.



Why Choose Strongarm Healthcare?



Customized Solutions

Every solution is customized to meet our clients' unique needs, ensuring seamless integration, exceptional ergonomics, and reliable performance.



Rapid Design and Prototyping

Our engineering team rapidly designs and prototypes custom solutions to meet demanding timelines without compromise.



Consignment Demo Product Program

We offer a generous 60-day consignment program, allowing you to evaluate our products in your own clinical or hospital environment at no cost or obligation.



Durability and Lower Cost of Ownership

Built for long-term performance, our products help reduce total cost of ownership, improve operational efficiency, and deliver lasting value—all backed by our industry-leading **limited lifetime warranty**.



Our products are proudly designed and built in Philadelphia, PA, USA

STRONGARM[®] + healthcare

Sales: sales@strongarmhealthcare.com

Support: support@strongarmhealthcare.com

 (800) 778-7901

 strongarmhealthcare.com



Strongarm Healthcare
425 Caredrean Drive
Horsham, PA 19044